

**INTENDED USE :**



This reagent kit is intended for “in vitro” quantitative determination of SGOT (AST) activity in serum/plasma.

**CLINICAL SIGNIFICANCE:**

Even though glutamate oxalate transaminase is widely distributed in various tissues of the body, it is a useful parameter in evaluating liver function. The elevated serum levels are found in case of hepatitis, obstructive jaundice, metastatic carcinoma, hepatic congestion and myocardial infarction or in kidney diseases

**PRINCIPLE:**

Kinetic determination of Aspartate aminotransferase (AST) based upon IFCC recommendations:

**REACTION:**

AST = Aspartate Aminotransferase

MDH = Malate dehydrogenase

**PACK SIZE & REAGENTS :**

**Reagent 1 :** Enzyme Reagent 1 x 10 ml 2 x 10 ml

**Reagent 2:** Substrate Reagent-ketoglutarate 5 x 10 ml 2 x 40 ml

**MATERIALS REQUIRED BUT NOT PROVIDED**

- Clean & Dry Glassware.
- Laboratory Glass Pipettes or Micro pipettes & Tips.
- Colormeter or Bio-Chemistry Analyzer.

**SAMPLES :**

Serum free of hemolysis. Heparin or EDTA plasma

**PREPARATION OF REAGENT & STABILITY:**

Mix 4 Volumes of Reagent 2 with One Volume of Reagent - 1.

Stability : 5 days at 20°C - 25°C, 4 weeks at 2°C - 8°C

Both the reagents are ready to use.

**GENERAL SYSTEM PARAMETERS**

Reaction type	: Kinetic reaction	Interval Time	: 60 Sec.
Wave length	: 340 nm	No. of readings	: 4
Cuvette Temperature	: 37°C	Factor	: 1746
Reagent Volume	: 1 ml	Zero setting	: Deionised water
Sample Volume	: 100	Light path	: 1cm
Delay Time	: 60 Sec	Units	: IU/L

**PROCEDURE:**

Working Reagent	1 ml
Sample	100 µl

Mix and after 60 second incubation, measure the change of Optical density per minute (AODI min.) during 4 minutes.

**CALCULATION :**

At 340nm with working reagent procedure for 1cm. Light path cuvette. Activity (IU/L)  $\Delta \text{Amin.} \times 1746$

**LINEARITY :**

Linear up to 300 IU/L

**NORMAL VALUE:**

Serum/Plasma < 40 IU/L

Normal values for infants are higher than adults.

Each Laboratory should establish its own expected range.

**QUALITY CONTROL:**

For accuracy it is necessary to run known controls with every assay.

**LIMITATION & PRECAUTIONS:**

1. Storage conditions as mentioned on the kit to be adhered.
2. Do not freeze or expose the reagents to higher temperature as it may affect the Performance of the kit
3. Before the assay bring all the reagents to room temperature.
4. Avoid contamination of the reagent during assay process.
5. Use clean glassware free from dust or debris.
6. Reagent: sample ratio as mentioned here above must be strictly observed as a change in it will adversely affect the factor.
7. Higher AST/GOT values may induce a false low result due to depletion of the substrate (total consumption of NADH before reading of the results). If an analyzer is used verify the presence of depletion factors on application.

**BIBLIOGRAPHY :**

1. Henderson A.R., Moss D.W., Enzymes, Tietz Fundamentals of clinical Chemistry, 5th Ed.
2. Vassault A., et al; Protocole de validation de techniques, Ann. Biol., clin. (1986) 44,686.